

Therapy - Minimize harm to patients



CHONGQING HAIFU MEDICAL TECHNOLOGY CO., LTD  
Tel: 86-23-6788 6788/6788 6199 Fax: 86-23-6788 6188  
Email: sales@haifu.cn Website: www.haifumedical.com  
Address: NO.1 Qingsong Road, Renhe, Yuebei District, Chongqing 401121, P.R. China  
Nov 2015



Haifu



Model JC  
Focused Ultrasound Tumor Therapeutic System

Model JC  
Focused Ultrasound Tumor Therapeutic System

#### Product Structures




Treatment Table  
• High-Frequency Generator  
• Integrated Transducer  
• 6-Dimension Motion Devices



Central Console  
• Ultrasound Monitoring Device  
• Therapeutic Control Part



Auxiliary Systems  
• Water Treatment System  
• Safety Protecting Device



#### Technology


- An ultrasound beam can propagate through living tissue harmlessly and be focused at a tiny focal region. The energy in the focal region is high enough to induce an immediate thermal toxicity (temperature above 56 °C) which will cause irreversible coagulative necrosis (a "lesion").

#### 3D Conformal Treatment

- From a point (single exposure) to a line, then a slice, then a volume that covers the entire tumor at any shape.
- Large-volume ablation in a single treatment
- Safe ablation of malignant tumors adjacent to major blood vessels

#### Powerful TPS Software

- 3D targeting module defines the boundary of tumor
- 3D planning module divides tumors into appropriate slices, records and analyzes coordinate information, forms a 3D therapeutic plan.
- Therapeutic module manages the treatment in conformity with treatment planning, monitors the tissue response and the safety of acoustic pathway, adjusts therapeutic parameters.



#### Precise Ablation Technology

##### Precise Boundary

- The treatment planning software enables conformal ablation of the whole tumor with no upper limit on the volume nor shape.
- The margin between treated and untreated tissue can be as narrow as 6 to 10 cells wide.

##### Precise Dosage

- Real-time imaging allows visual feedback during treatment process
- Immediate image after each exposure can be compared with the previous
- An operator can adjust the dose anytime to suit the individual needs
- Integrated dose data will be recorded for future analysis and effect evaluation

##### Precise Control

- With 11 mm accumulative error, the accurate movement of 6-dimensional motion system can ablate tumors adjacent to major vessels and nerves safely
- Color Doppler Ultrasound provides clear real-time monitoring during the whole treatment procedure



Indications


Solid tumors of soft tissues, such as liver tumor, uterine fibroids, breast cancer, kidney tumor, pancreas tumor, bone tumor.

Contraindications


- Tumors in air-containing viscera, such as lung, stomach, and bowel
- Mediastinal tumors
- Spinal tumors

Clinical Advantages


- Noninvasive treatment to preserve organs and structures, with no blood transfusion and no radiation
- Conformal and precise ablation
- One time treatment, no limitation by tumor size and shape
- Real-time ultrasound-guided therapy with digital quantitative analysis
- Activation of immune system



Lian Zhang, High-intensity focused ultrasound (HIFU) effective and safe therapy for hepatocellular carcinoma adjacent to major hepatic veins. Eur Radiol (2009) 19:437-445



Peng Wei, High-intensity focused ultrasound ablation of breast cancer. Expert Rev. Anticancer Ther. 7(6), (2007), 823-831.



Wenzhi Chen, Primary Bone Malignancy Effective Treatment with High-Intensity Focused Ultrasound Ablation, Radiology: Volume 235: Number 3—June 2010

Liver cancer


A 51-year-old man with HCC. (a) Before HIFU, CT image shows a large tumor adjacent major hepatic veins (b) Enhanced MR images shows no enhancement in treated region and no damage to hepatic vein 2 weeks after HIFU treatment.

Breast cancer


Contrast-enhanced MRI shows complete necrosis of HIFU-treated breast cancer. (a) Before treatment, the breast lesion is circumscribed, and contrast uptake is detected in the tumor (arrow) and normal tissue surrounding the tumor. (b) 7 days after HIFU, an absence of contrast uptake in ablated volume (arrows), including tumor and a margin of treated normal breast tissue about 1.5–2.0 cm around the cancer.

Bone tumour

Enhanced MR images obtained in 18-year-old woman who underwent HIFU ablation for tibia osteosarcoma. (a) Image before HIFU shows single hypervascular lesion (arrow) in the left tibia. Images (b) 2 weeks and (c) 36 months after HIFU show no evidence of enhancement in treated region (arrow).



Gianluigi Orgera, High Intensity Focused Ultrasound Ablation of Pancreatic Neuroendocrine Tumours: Report of Two Cases. Cardiovasc Intervent Radiol 2011 Apr;34(2):419-423



Pancreatic cancer

A 43-year-old woman with pancreatic tumor. (a) MRI before HIFU shows a 2.4-cm lesion (arrowheads). (b) MDCT performed 24h after HIFU shows ablation of entire lesion (arrowheads). (c) MDCT scan 9 months after treatment shows the lesion appears slightly diminished in dimensions but widely hypodense and without residual enhancing region (arrowheads).

Uterine fibroids

Enhanced MRI images in a 37-year-old woman with uterine fibroids. Before treatment, the fibroid had abundant and active blood supply. Fourteen days after treatment, the fibroid showed obviously shrinkage and lost activity. MRI examination 4.5 months and 10 months after HIFU showed obvious shrinkage of treated fibroids.

End-users

Well-known hospitals in the world are using Hailufu's system:



Over 100 centers and 50,000 cases worldwide (Feb. 2013)

European Institute of Oncology, Milan, Italy

- European Institute of Oncology, Milan, Italy
- The John Radcliffe Hospital, Oxford, UK
- University Clinics of Bonn, Bonn, Germany
- Hospital Mutua de Terrassa, Barcelona, Spain
- Saint Marina Hospital, Plovdiv, Bulgaria
- Medical Center of Central Bank of Russian Federation, Moscow, Russia
- Queen Mary Hospital of the University of Hong Kong, China
- King Fahad Medical City, Riyadh, Saudi Arabia
- Incheon Christian Hospital, Incheon, Korea
- CIMES (Centro de Investigaciones Médico Quirúrgicas), Havana, Cuba
- 301 PLA General Hospital, Beijing, China
- etc.

Total Solution

Professional Equipment

With complete intellectual rights and CE-mark, it is the world first equipment clinically applied in tumor treatment, a result of over 20 years experience and optimization.



Experienced Specialists

A team of experienced medical and engineering specialists will provide integrated training and service to enable the end-users independent operation of the equipment.



Customized Solution

Suitable clinical protocols, operation & management, advices and research cooperation proposal will be tailored for each end-user.



Qualification



|                          |  |  |
|--------------------------|--|--|
| Main Parameters          | Acoustic focusing efficiency             | 28000  |
|                          | Focal region                             | 1.1mm×1.1mm×3.3mm  |
|                          | Max acoustic intensity                   | 25000W/cm <sup>2</sup>   |
|                          | Max output acoustic power                | 400W   |
|                          | Side lobe                                | <10dB  |
|                          | Maximum range of transducer movement     | X=120mm,Y=120mm,Z=180mm  |
|                          | Movement control accuracy                | ±0.1mm   |
|                          | Accumulated tolerance in linear movement | ±1mm   |
|                          | Therapeutic frequency                    | 0.8-2.4MHz   |
|                          | Probe vertical movement range            | 0-100mm  |
|                          | Probe rotating angle range               | ±90°   |
|                          | Dissolved oxygen                         | ≤3ppm  |
|                          | Electrical Power                         | 12.5KVA  |
| Installation Environment | Room requirement                         | Area : ≥40m <sup>2</sup> , Width : ≥5m   |
|                          | Power requirement                        | Three-phase five-wire power cable with ground wire which in conformity with local laws |
|                          | Water requirement                        | Flow: ≥1m <sup>3</sup> /h, Pressure:0.1-0.5MPa   |

